

IMAGING – ADDITIONAL NOTES FOR LSE/HUGHES GUN

The current for each aperture is different than for UMER. Measure each using current monitor.

LSE Station:

The phosphor screen, camera and software is different from the IMPERX we commonly use on UMER:

- A. We use fast phosphor screen with about 3-4ns response/decay
- B. We use a microchannel plate (MCP) gated intensified CCD (ICCD) camera (PIMAX2) which is capable of producing 16bit images (see: www.piacton.com for description of PIMAX2 camera)

We will examine two gun aperture settings

- 1) 25mA aperture (equivalent to 22mA aperture on UMER)
- 2) Five beamlet aperture

Pay attention to changes in beam current trace as a function of BIAS voltage: use: 30, 40, 45 V bias.

Perform beam imaging experiments

- 1) with camera gate > 100 ns
- 2) one case with short camera gate < 10 ns to get a time-resolved image

Imaging Procedure: Instructor will demonstrate camera turn on, software enable, gate mode enable

Student data acquisition:

- 1) SET Gates per exposure to 1
- 2) SET Gain to 1 (0-256 nonlinear)
- 3) SET Gate width (100ns nominal)
- 4) SET Gate delay (318ns nominal)
- 5) Check oscilloscope to see if gate overlaps Current Pulse
If NOT adjust Gate width and/or Gate delay
- 6) Enable FOCUS MODE
- 7) Adjust camera lens focus ring to get well defined beam image
- 8) Adjust Gain, Gates per exposure to get about 5000 counts in center of image
- 9) STOP FOCUS MODE
- 10) ENABLE ACQUISITION MODE
- 11) SAVE image as 8 bit TIFF in directory USPAS_STUDENT_IMAGES and
DOWNLOAD TO MEMORY STICK for processing